

Titanium Screw Guide

IMPLANT SYSTEM	Prosthetic Screw Connection	Prosthetic Screw Thread	Mfg. Recommended Torque
3i Certain Internal	0.048" (1.22mm) Hex	M1.6 x 0.35	20 Ncm
3i Certain External	Square	M2 x 0.40	32-35 Ncm
Ankylos C/X	1mm Hex	M1.8 x 0.35	15 Ncm
Astra OsseoSpeed			
X-Small	0.050" (1.27mm) Hex	M1.4 x 0.30	15 Ncm
Small	0.050" (1.27mm) Hex	M1.6 x 0.35	20 Ncm
Large	0.050" (1.27mm) Hex	M2 x 0.40	25 Ncm
BioHorizon Internal/External	0.050" (1.27mm) Hex		30 Ncm
Brånemark			
Brånemark System NP	Unigrip	M1.6 x 0.35	35 Ncm
Brånemark System RP	Unigrip	M2 x 0.40	35 Ncm
Brånemark System WP	Unigrip	M2.5 x 0.45	35 Ncm
Camlog ScrewLine			
3.3, 3.8, 4.3mm	0.050" (1.27mm) Hex	M1.6 x 0.35	20 Ncm
5.0, 6.0mm	0.050" (1.27mm) Hex	M2 x 0.40	20 Ncm
Endopore (Sybron)	0.050" (1.27mm) Hex		30 Ncm
Inclusive Tapered	0.050" (1.27mm) Hex	1-72 UNF	35 Ncm
Keystone Prima Connex	Square	M1.8 x 0.35	30 Ncm
Keystone Restore			30 Ncm
SD Restore	0.048" (1.22mm) Hex, Square	1-72 UNF	30 Ncm
RD Restore	0.048" (1.22mm) Hex, Square	M2 x 0.40	30 Ncm
WD Restore	0.048" (1.22mm) Hex, Square	M2.5 x 0.45	30 Ncm
Keystone Renova	0.048" (1.22mm) Hex	M1.8 x 0.35	30 Ncm
MIS Seven	0.050" (1.27mm) Hex	1-72 UNF	30 Ncm
Neoss	Neoss Proprietary	M2 x 0.40	32 Ncm
NobelActive			
NobelActive 3.0	Unigrip	M1.4 x 0.30	15 Ncm
NobelActive NP	Unigrip	M1.6 x 0.35	35 Ncm
NobelActive RP/WP	Unigrip	M2 x 0.40	35 Ncm
NobelReplace			
NobelReplace NP	Unigrip	M1.8 x 0.35	35 Ncm
NobelReplace RP	Unigrip	M2 x 0.40	35 Ncm
OCO Biomedical	0.050" (1.27mm) Hex	1-72 UNF	30 Ncm
Straumann Tissue-Level			
Narrow Neck	SCS	M1.8 x 0.35	35 Ncm
Regular Neck	SCS	M2 x 0.40	35 Ncm
Wide Neck	SCS	M2 x 0.40	35 Ncm
Straumann Bone-Level			
Narrow Cross-fit	SCS	M1.6 x 0.35	35 Ncm
Regular Cross-fit	SCS	M1.6 x 0.35	35 Ncm
Sybron Pro			
Sybron Pro (Internal Octa)	0.050" (1.27mm) Hex	M2 x 0.40	35 Ncm
Sybron Pro (Internal Hex)	0.050" (1.27mm) Hex	M2 x 0.40	35 Ncm
Zimmer Spline	0.050" (1.27mm) Hex	1-72 UNF	30 Ncm
Zimmer Screw-Vent	0.050" (1.27mm) Hex	1-72 UNF	30 Ncm